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4-Step Procedure for Laboratory Cylinder Acceptance


(Excluding the dissolved acetylene cylinders and cryogenic insulated cylinders)

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
Inspect the appearance of the cylinder (the gas inside matches the color of the cylinder, the cylinder outlet is sealed with a plastic film or plug, no deformation, crack, abrasion, corrosion of the cylinder and oxygen cylinders are free from grease) and verify the integrity of safety accessories (shock-absorbing rubber ring, cylinder protective cover, valve).

Scan the electronic identification symbol on the cylinder shoulder (such as a QR code) using a mobile phone to view basic information about the cylinder, including cylinder status and specific expiration date of the inspection.

Check the manufacturer's certificate of compliance (including gas name, purity, volume, filling date, expiration date, manufacturer's name, etc.) and labeling.

Check the steel mark on the shoulder of the cylinder (including the manufacture date, the designed service life, the name or the chemical formula of the filled gas, etc.).

** If there are any problems or concerns in the above four steps, please contact the cylinder supplier or feedback to LSMD in time. Make sure there is no any problem before accepting the cylinder.